

NWL lichen ID kit, without reagents, July 2023

P and K reagents can be purchased in dry form at NWL meetings, workshops, and trainings

Number	Item
1	Forceps, straight, ultra fine, stainless steel
1	Teasing needle (dissecting needle, steel)
1	Beaker, 25 ml, Pyrex
1	Spotting plate, 12-well, ceramic
12	Double edged razors, Personna
4	Bottles, amber glass, with droppers and lids (for distilled water, ethanol, acetone, and liquid of your choice, e.g. IKI; reagents not included)
2	Bottles, 15 ml, polyethylene, with dropper cap (for K and C)
11	Microscope slides, glass
30	Cover slips, glass, for micro slides
16	Capillary tubes, glass, nonheparinized, 1.1 X 75 mm
1	Protocol for applying P reagent (p-phenylenediamine)
2	Chemical safety data sheets for P, (KOH), and bleach; OSHA HCS
1	Polyethylene container for kit with lid
